


White Glue

Light-Cure

Orthodontic Adhesive


Elegant Aesthetics

White Glue keeps clear white by clinical treatment which is similar to tooth color. It meets the aesthetic demands of people, which can be used to bond orthodontic metal/ ceramic/ sapphire brackets, buccal tubes and molar bands.



Strong Bonding Strength

Bonding Strength of adhesive can reach 20MPa (min. Clinical Effective Bonding Strength is 7MPa), which assures the bonding stability of brackets.



High-efficient Treatment

White Glue can save the chair time greatly to have better treatment result.

- ▶ Dentists have enough operation time to locate brackets and remove excess adhesive under ambient light.
- ▶ The light curing time can be shortened to 5 seconds.



Excellent Experience

- ▶ **Etchant Gel:** Gel-type without overflowing, moderate viscosity
- ▶ **Primer:** Moisture-tolerant, good fluidity, strong permeability, without irritating odor
- ▶ **Adhesive:** Fluent extrusion, clear and smooth

White Glue | Light-Cure Orthodontic Adhesive Kit

- ✂ ×2 Adhesive 4.0g/piece
- ✂ ×1 Etchant 3.0ml/piece
- 🍷 ×1 Primer 3.0ml/piece
- 🪥 ×20 Brush
- 🪡 ×2 Syringe



Hangzhou Xingchen 3B Dental Instrument & Material Co.,Ltd.



Add: Building 4, No.16 Xingyuan Road, DaYuan Town, FuYang District, 311421, Hangzhou, China

Tel: +86-571-87655303

E-mail: hz3b.n1@hz3b.com

Website: www.ortho3b.com